

## Powered stair climbing evacuation & rescue chair...



### What's an Escape-Chair® VOLT?

An Escape-Chair® VOLT is the **POWERED** solution to evacuate & rescue people to safety in a seated position up and/or down a flight of stairs. A simple procedure can make the Escape-Chair® VOLT fully functional within a few seconds.

### How does the Escape-Chair® VOLT work?

The Escape-Chair® VOLT is operated by a self-explanatory control panel with display on the operating handle. After unfolding the Escape-Chair® VOLT it is immediately ready for use. The single operator must push the operating handle down in the direction of the stairs.

### Advantages of the Escape-Chair® VOLT:

- Powered support makes going up the stairs much easier for operator.
- Ideal for quick evacuation (3 speed descent).
- Lithium cell enough for up/down 100 floors.
- Single operator chair.
- Ready for use within seconds.
- Easy to use and operate (wireless).
- Equipped with arm- and footrest.
- 3 height adjustable operating handle.
- Flame & heat resistant seat, headrest and dust cover.
- Battery indicator on display control panel.
- LED light for dark environment evacuation.
- CE certified.

### Our advice:

Make sure you have an Escape-Chair® VOLT available in your organization so that you're prepared for any evacuation and/or emergency, for up and/or down evacuation.

"Failure to prepare for an emergency has its consequences."





Specifications	Escape-Chair® VOLT
Article number:	10050
Dimensions folded up (height x width x depth):	115 x 50 x 27 cm
Maximum load bearing:	180 kg
Maximum stair gradient:	40°
Weight:	±30 kg
Using time when full charged:	2 hours
Motor:	36V, 300W, no noise
Repeated charging times (on 220Volt):	500
Battery life:	3 - 5 years
Extra batterypack:	Optional
Number of speeds:	3
Ready for use within seconds:	✓
Foot rest:	✓
Arm rests:	✓
LED lightning for dark environment:	✓
3 height adjustable operating handle:	✓
Flame & heat resistant seat, headrest & dust cover:	✓
Carrying handles front & back:	✓
CE certified:	✓
Free demonstration and advise provided:	✓
Delivery incl. manual, identity sign & wall hooks:	✓
Warranty on non-wearing parts:	1 year

Interim changes to technical specifications and printing errors reserved.

