

Escape-Sheet® Economy Line

Evacuation directly ... from a bed

What`s an Escape-Sheet®?

People who aren't able to leave their bed independently need to be evacuated using evacuation equipment. One of these aids is the Escape-Sheet®. A mattress with an Escape-Sheet® underneath can be transformed very quickly into an evacuation mattress.

How does an Escape-Sheet® work?

In case of emergency the strong safety belts are taken from underneath the mattress. The belts are strapped over the patient and bedding so the patient is safely secured to the mattress. The Escape-Sheet® can then be easily pulled to the floor and taken to a place of safety.

Advantages Escape-Sheet®:


- Plastic buckle pockets increase the life expectancy of the product (allocated on the side of the sheet).
- Handgrips to assist with sideways manoeuvrability.
- Suitable for each mattress height.
- TÜV tested.
- CE certified.

Our advice:

Make sure every bed in your Organization is equipped with an Escape-Sheet®, so you're always prepared for optimal and rapid evacuation.

A LABEL
QUALITY
EVACUATION
PRODUCT





Model Economy Line



Washing label



Suitable for all mattress heights



Buckles in pockets

Specifications	Escape-Sheet® Economy Line	CE
Article number:	21270	
Dimensions (length x width):	± 200 x 90 cm	
Material:	67% Polyester / 33% Cotton	
Maximum payload:	200 kg	
Weight:	ca. 800 gram	
Demonstration and advise provided:	✓	
Ready within seconds:	✓	
Immediate evacuation directly from the bed:	✓	
Industrial wash 60° C:	✓	
For all mattress thicknesses:	✓	
Buckles in pockets:	✓	
Warranty:	2 years	



Interim changes to technical specifications and printing errors reserved.
 Produced by Escape Mobility Production & Innovation: TÜV approved organization.
 Health guaranteed through Öko-tex Standard 100

