

Escape-Chair® VOLT



Powered stair climbing evacuation & rescue chair...

What's an Escape-Chair® VOLT?

An Escape-Chair® VOLT is the **POWERED** solution to evacuate & rescue people to safety in a seated position up and/or down a flight of stairs. A simple procedure can make the Escape-Chair® VOLT fully functional within a few seconds.

How does the Escape-Chair® VOLT work?

The Escape-Chair® VOLT is operated by a self-explanatory control panel with display on the operating handle. After unfolding the Escape-Chair® VOLT it is immediately ready for use. The single operator must push the operating handle down in the direction of the stairs.

Advantages of the Escape-Chair® VOLT:

- Powered support makes going up the stairs much easier for operator.
- Ideal for quick evacuation (4 speed descent).
- Lithium cell enough for up/down 40 floors.
- Single operator chair.
- Ready for use within seconds.
- Easy to use and operate (wireless).
- Equipped with arm- and footrest.
- 4 height adjustable operating handle.
- Battery indicator on display control panel.
- LED light for dark environment evacuation.
- CE certified.

Our advice:

Make sure you have an Escape- Chair® VOLT available in your organization so that you're prepared for any evacuation and/or emergency, for up and/or down evacuation.

"Failure to prepare for an emergency has its consequences."





| Specifications | Escape-Chair® VOLT |
|--|-----------------------------------|
| Article number: | 10050 |
| Dimensions folded up (height x width x depth): | 109 x 63 x 29 cm |
| Maximum load bearing: | 180 kg |
| Maximum stair gradient: | 40° |
| Weight: | ±37,5 kg |
| Using time when full charged: | 2 hours (± 40 floors) |
| Type of tracks: | Serrated tracks for superior grip |
| Unlocking function tracks: | ✓ |
| Battery life: | 2 – 3 years |
| Extra batterypack: | Optional |
| Number of speeds: | 4 |
| Ready for use within seconds: | ✓ |
| Foot rest: | ✓ |
| Arm rests: | ✓ |
| LED lightning for dark environment: | ✓ |
| 4 height adjustable operating handle: | ✓ |
| Flame & heat resistant dust cover: | Optional |
| Carrying handles front & back: | ✓ |
| 10G crash test bracket for Ambulances / vehicles | ✓ |
| CE certified: | ✓ |
| Wall mounting bracket included | ✓ |
| Delivery incl. wall bracket: | ✓ |
| Warranty on non-wearing parts: | 1 year (Battery 6 months) |

Interim changes to technical specifications and printing errors reserved.

