



Escape Chair VOLT



Powered stair climbing evacuation & rescue chair...

What's an Escape-Chair VOLT?

An **Escape-Chair VOLT** is the **POWERED** solution to evacuate & rescue people to safety in a seated position up and/or down a flight of stairs.

A simple procedure can make the **Escape-Chair VOLT** fully functional within a few seconds.

How does the Escape-Chair VOLT work?

The **Escape-Chair VOLT** is operated by a self-explanatory control panel with display on the operating handle. After unfolding the **Escape-Chair VOLT** it is immediately ready for use. The single operator must push the operating handle down in the direction of the stairs.

Advantages Escape-Chair VOLT:

- Powered support makes going up the stairs much easier for operator.
- Ideal for quick evacuation (4 speed descent).
- Lithium cell enough for up/down 40 floors.
- Single operator chair.
- Ready for use within seconds.
- Easy to use and operate (wireless).
- Equipped with arm- and footrest.
- 3 height adjustable operating handle.
- Flame & heat resistant seat, headrest and dust cover.
- Battery indicator on display control panel.
- LED light for dark environment evacuation.
- Big wheels for a more comfortable ride
- CE & ISO certified.

Specifications	Escape-Chair VOLT
Article number:	10050
Dimensions folded up (height x width x depth):	109 x 63 x 29 cm
Maximum load bearing:	254 kg
Maximum stair gradient:	45°
Weight:	±37 kg
Using time when full charged:	2 hours (±40 floors)
Battery life:	2-3 years
Extra battery pack:	Optional
Number of speeds:	4
Ready for use within seconds:	✓
Anti-blocking system on motor:	✓
Extra-large wheels for comfort and easy handling:	✓
Foot rest & arm rests:	✓
LED lightning for dark environment:	✓
Battery indicator on display:	✓
4 height adjustable operating handle:	✓
Carrying handles front & back:	✓
CE & ISO certified:	✓
User manual, identity sign & wall hooks:	✓
Warranty on non-wearing parts:	10 years

Our advice:

Make sure you have an **Escape-Chair VOLT** available in your organization so that you're prepared for any evacuation and/or emergency, for up and/or down evacuation.

"Failure to prepare for an emergency has its consequences."



E: info@escape-mobility.com

Escape Mobility Company Ltd.
85 Great Portland Street, 1st floor,
London, W1W 7LT

T: +44 (0)1908 229 028

W: www.escape-mobility.com

Escape Mobility International BV.
Jan Campertstraat 7 – Unit 1.14,
6416 SG, Heerlen.

T: +31 (0)45 528 0670

...your way to safety

