



Escape-Mattress Compact



Best choice for healthcare environment and public buildings

What's an Escape-Mattress?

Escape-Mattress models are designed for individuals who would not be able to comfortably get into an evacuation chair, or where evacuation chairs may struggle (for example when evacuating an individual down a spiral staircase). Evacuation mattresses are ideal for medical centres and hospitals, as they cater for complex healthcare needs. E.g. bed-bound and unconscious patients. **The Escape-Mattress** is an "A-label quality" product and produced in-house.

Our mattresses are CE, TÜV NORD and ISO certified for up to 150kg payload.

We guarantee the highest quality standards.

The **Escape-Mattress Compact** is a complete **Escape-Mattress** but its smaller dimensions allows it to be used more flexibly. The **Escape-Mattress Compact** is rolled up and stored in the storage bag that is included. This **Escape-Mattress Compact** is more suitable, so that it can be taken to different situations where an evacuation aid is required.

The **Escape-Mattress** is made of fire resistant material that is resistant to blood, urine, acids and alkaline as well.

CE, TÜV NORD and ISO certified for up to 150kg payload.

There is no better and more durable quality available on the market than any of our **Escape-Mattress** range.

Advantages of the Escape-Mattress Compact:

- Stored double-folded (space saving)
- 4 carrying handles on each side

Protection of evacuee/patient:

- 2 body flaps – 1 per side
- Foot bag
- Adjustable headrest

Specifications Escape-Mattress Compact	
Article number:	21030
Ready to use within a few seconds:	✓
Heat and flame retardant:	✓
Dimensions (length x width):	200 x 68 cm
Mattress thickness:	6 cm
Including storage bag or case	✓
Maximum carrying capacity:	150 kg
Space-saving:	Can be rolled up
Weight:	7.5 kg incl. bag
Adjustable headrest with immobilisation band:	✓
Bag dimensions:	80 x 70 x 21 cm
Blood, urine, & acid resistant:	✓
5-language manual:	✓
Photoluminescent identity sign:	✓
CE, TÜV NORD and ISO certified:	✓
Warranty:	2 years

Usage:

The **Escape-Mattress Compact** is mainly used in healthcare environments and (public) buildings.



E: info@escape-mobility.com

Escape Mobility Company Ltd.
85 Great Portland Street, 1st floor,
London, W1W 7LT

T: +44 (0)1908 229 028

W: www.escape-mobility.com

Escape Mobility International BV.
Jan Campertstraat 7 – Unit 1.14,
6416 SG, Heerlen.

T: +31 (0)45 528 0670

...your way to safety

