



Escape-Mattress Economy Line

The **Escape-Mattress Economy Line** can be triple folded so that it is easy to move and store. Once unfolded, this evacuation mattress is immediately ready to use. This mattress can be stored both hygienically and dust-free, by using the storage bag or the Escape-Locker.



Usage:

The **Escape-Mattress Economy Line** is mainly used in (public) buildings and in Care Homes.

Best choice for healthcare environment and public buildings

What's an Escape-Mattress?

Escape-Mattress models are designed for individuals who would not be able to comfortably get into an evacuation chair, or where evacuation chairs may struggle (for example when evacuating an individual down a spiral staircase). Evacuation mattresses are ideal for medical centres and hospitals, as they cater for complex healthcare needs. E.g. bed-bound and unconscious patients. **The Escape-Mattress** is an "A-label quality" product and produced in-house.

Our mattresses are CE, TÜV NORD and ISO certified for up to 150kg payload.

We guarantee the highest quality standards.

Advantages Escape-Mattress Economy Line:

- Ready to use within seconds
- Space-saving: triple folding
- Can be stored in Escape-Locker
- Storage bag included
- 2 safety straps
- Foot bag with zipper
- Adjustable head rest

Specifications Escape-Mattress Economy Line

Article number:	21040
Ready to use within a few seconds:	✓
Dimensions (length x width):	200 x 65 cm
Material:	100% polyester with PVC coating
Mattress thickness:	5 cm
Including storage bag or case	✓
Maximum carrying capacity:	150 kg
Space-saving:	Can be wall mounted
Weight:	4.8 kg
Adjustable headrest with immobilisation band:	✓
Warranty:	2 years
5-language manual:	✓
CE, TÜV NORD and ISO certified:	✓



E: info@escape-mobility.com

Escape Mobility Company Ltd.
85 Great Portland Street, 1st floor,
London, W1W 7LT

T: +44 (0)1908 229 028

W: www.escape-mobility.com

Escape Mobility International BV.
Jan Campertstraat 7 – Unit 1.14,
6416 SG, Heerlen.

T: +31 (0)45 528 0670

...your way to safety

