



Escape-Mattress Junior



Best choice for maternity departments, playgroups and daycare centres.

What's an Escape-Mattress?

Escape-Mattress models are designed for individuals who would not be able to comfortably get into an evacuation chair, or where evacuation chairs may struggle (for example when evacuating an individual down a spiral staircase). Evacuation mattresses are ideal for medical centres and hospitals, as they cater for complex healthcare needs. E.g. bed-bound and unconscious patients. **The Escape-Mattress** is an "A-label quality" product and produced in-house.

Our mattresses are CE, TÜV NORD certified & ISO tested for up to 150kg payload.

We guarantee the highest quality standards.

Usage:

The **Escape-Mattress Junior** is excellent for the safe and quick evacuation of maternity departments, playgroups and daycare centers.

The **Escape-Mattress Junior** is designed for evacuating 2 toddlers (up to 3 years) at the same time in lying position. The **Escape-Mattress Junior** can be used by a single operator only. This means that the toddlers do not have to be carried down the stairs in the event of an emergency evacuation.

The **Escape-Mattress Junior** is stored in the supplied storage bag. CE, TÜV NORD certified & ISO tested for up to 150kg payload.

There is no better and more durable quality available on the market than any of our **Escape-Mattress** range.

Advantages of the Escape-Mattress Junior:

- Ready to use within seconds by just 1 operator only.
- Evacuating 2 toddlers (up to 3 years) at the same.
- Space saving (wall mounting).
- Ideal for maternity ward, day care and kindergarten.
- Made of high quality & fire resistant materials.
- 1 large and 1 small tow belt at both head and foot ends.
- Blood-, urine-, acid and alkali resistant.
- 2 adjustable head rest.
- 2 body bags.
- Comes with storage bag.
- CE, TÜV NORD certified & ISO tested.

Specifications Escape-Mattress Junior	
Article number:	21070
Dimensions (length x width x thick):	± 120 x 68 x 8 cm
Maximum payload:	150 kg
Age evacu�e(e)	Up to 3 years
Weight:	± 6 kg incl. cover
Heat and flame retardant:	✓
Ready for use within seconds:	✓
Speed of descent to be determined by the user:	✓
Blood, urine & acid resistant:	✓
Equipped with headrest & body-covers:	✓
Including storage bag:	✓
Photoluminescent identity sign:	✓
Wall hooks:	✓
5-language manual:	✓
CE, TÜV NORD certified & ISO tested:	✓
Warranty:	2 years



E: info@escape-mobility.com

Escape Mobility Company Ltd.
85 Great Portland Street, 1st floor,
London, W1W 7LT

T: +44 (0)1908 229 028

W: www.escape-mobility.com

Escape Mobility International BV.
Jan Campertstraat 7 – Unit 1.14,
6416 SG, Heerlen.

T: +31 (0)45 528 0670

...your way to safety

