



Escape-Mattress Solo



Best choice for healthcare environment and public buildings

What's an Escape-Mattress?

Escape-Mattress models are designed for individuals who would not be able to comfortably get into an evacuation chair, or where evacuation chairs may struggle (for example when evacuating an individual down a spiral staircase). Evacuation mattresses are ideal for medical centres and hospitals, as they cater for complex healthcare needs. E.g. bed-bound and unconscious patients. The **Escape-Mattress** is an "A-label quality" product and produced in-house.

Our mattresses are CE, TÜV NORD & ISO certified up to 150kg payload.

We guarantee the highest quality standards.

The **Escape-Mattress Solo** is fastened to the wall. By opening the zip at the underside of the cover, the **Escape-Mattress Solo** slides out and is ready to use. Suitable for places where there is limited space, such as operating theatres.

The **Escape-Mattress** is made of fire resistant material that is resistant to blood, urine, acids and alkaline as well.

Our mattresses are CE, TÜV NORD & ISO certified up to 150kg payload. There is no better and more durable quality available on the market than any of our **Escape-Mattress** range.

Advantages Escape-Mattress Solo:

- Supplied with a storage cover
- Can be suspended from the wall
- 4 handles at each side

Protection of evacuee/patient:

- 2 body flaps – 1 per side
- Foot bag
- Adjustable headrest

Specifications Escape-Mattress Solo	
Article number:	21020
Ready to use within a few seconds:	✓
Heat and flame retardant:	✓
Dimensions (length x width):	200 x 68 x 8cm
Also available in a 55cm wide version:	✓
Mattress thickness:	8 cm
Including storage bag:	✓
Maximum carrying capacity:	150 kg
Space-saving:	Wall mounted
Weight:	9.5 kg incl. cover
Adjustable headrest with immobilisation band:	✓
Blood, urine & acid resistant:	✓
Photoluminescent identity sign:	✓
Wall hooks:	✓
5-language manual:	✓
CE, TÜV NORD & ISO certified:	✓
Warranty:	2 years

Usage:

The **Escape-Mattress Solo** is ideal for healthcare environments like operating theatres.



E: info@escape-mobility.com

Escape Mobility Company Ltd.
85 Great Portland Street, 1st floor,
London, W1W 7LT

T: +44 (0)1908 229 028

W: www.escape-mobility.com

Escape Mobility International BV.
Jan Campertstraat 7 – Unit 1.14,
6416 SG, Heerlen.

T: +31 (0)45 528 0670

...your way to safety

