



Escape-Sheet Custom Made



The **Escape-Sheet** is positioned permanently between the mattress and the bed frame. In the event of an emergency, the safety straps of the **Escape-Sheet** are secured over the mattress and the patient.

The patient is immobilised safely for transporting down the stairs. The **Escape-Sheet** is used frequently in hospitals and nursing homes.

The **Escape-Sheet** is widely used by hospitals and healthcare facilities. Assists in moving a patient to a place of safety, horizontally as well as vertically. Our sheets are water repellent, anti-static and also available in heat and flame resistant.

Advantages of the Escape Sheet Custom Made:

- Straps and handles highly visible in bright green colour.
- Extra-long towing straps.
- Plastic buckle pockets increase the life expectancy of the product
- Handgrips to assist with sideways manoeuvrability.
- Static-Control™ (anti-static polyester).
- Adjusting corner straps (suitable for each mattress height).
- Hydro-tech Nano technology (water-repellent)

Evacuation directly.... from a bed

What's an Escape Sheet?

People who aren't able to leave their bed independently need to be evacuated using evacuation equipment. One of these aids is the **Escape-Sheet**. A mattress with an **Escape-Sheet** underneath can be transformed very quickly into an evacuation mattress.

How does an Escape Sheet work?

In case of emergency the strong safety belts are taken from underneath the mattress. The belts are strapped over the patient and bedding so the patient is safely secured to the mattress. The **Escape-Sheet** can then be easily pulled to the floor and taken to a place of safety.

Our sheets are CE certified & TÜV certified up to 250kg. We guarantee the highest quality standards.

Specifications Escape Sheet Custom Made	
Article number:	21260
Dimensions (length x width):	On demand
Material:	99% Polyester, 1% Static-Control™
Maximum payload:	250 kg
Weight:	On demand
Corner strap adjusters:	✓
Industrial wash 75° C:	✓
Immediate evacuation directly from the bed:	✓
Ready to use within seconds	✓
Hydro-tech nano technology (water-repellent):	✓
Free demonstration and advise:	✓
For all mattress heights	✓
Body covers:	On demand
Extra cover at foot end:	On demand
Buckles in pockets:	On demand
Warranty:	2 years

Our advice:

Make sure every bed in your organization is equipped with an **Escape-Sheet**, so you're always prepared for optimal and rapid evacuation.



E: info@escape-mobility.com

Escape Mobility Company Ltd.
85 Great Portland Street, 1st floor,
London, W1W 7LT

T: +44 (0)1908 229 028

W: www.escape-mobility.com

Escape Mobility International BV.
Jan Campertstraat 7 – Unit 1.14,
6416 SG, Heerlen.

T: +31 (0)45 528 0670

